Stethoskope HEINE GAMMA 2.1®, 2.2®, 2.3®, GAMMA C®

CE 93/42/EWG/CEE/EEC

HEINE

HEINE Optotechnik GmbH & Co. KG
Kientalstr. 7 · 82211 Herrsching · Germany
Tel. +49(0)8152/38-0 · Fax +49(0)8152/38-202
E-Mail: info@heineopto.de · www.heine.com
med 3101 1/11.05

DEUTSCH

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Stethoskope HEINE GAMMA 2.1®, 2.2®, 2.3®, GAMMA C®

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Instructions for use

Stethoscopes HEINE GAMMA 2.1®, 2.2®, 2.3®, GAMMA C®

Read these instructions carefully before using the HEINE stethoscope and keep them in a safe place for future reference.

Setting-up

Position the chestpiece without applying too much pressure. It is not necessary to cover the whole surface of the membrane. The combi-chestpiece lets you choose between a membrane and an open cone. To select, hold the tube in one hand and rotate the chestpiece with the other hand 180° until it clicks into place.

Cleaning and maintenance

The ear tube and plastic tubing can be wiped clean with alcohol or mild soapy water.

The ear olives can be unscrewed and washed in soapy water.

The membrane can be removed after carefully unscrewing the retaining ring and washed in soapy water. A replacement membrane is supplied.

All components can be sterilized in ethylene oxide. Do not use thermal sterilization or soaking.

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Warranty

We guarantee this device for five years from date of purchase. We guarantee the proper function of the unit provided it is used as intended by the manufacturer and according to the instructions for use. Any faults or defects which occur during the guarantee period will be repaired free of charge provided they are caused by faults in material, design or workmanship. In the case of a complaint of a defect in the product, the buyer must prove that any such defect was present at the time the product was delivered. This guarantee excludes any defects which occur because of wear and tear, incorrect use, the use of non-original HEINE parts. Also excluded are repairs or modifications made by persons not authorised by HEINE or cases where the customer does not follow the instructions supplied with the product. Further claims, in particular claims for damage not directly related to the HEINE product, are excluded.

For repairs, please contact your supplier, who will return the goods to us or our authorized agent.

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